

NEW!!

Inner Glass Door with Access Windows for 30-Liter Incubator (Optional)



For SCA/SMA-30 Series: P/N GAW-30DH



For APC/APM-30 Series: P/N GAW-30AQ

Advantages of Inner Glass Door with Access Windows:

- Removable gaskets – easier cleaning
- Interchangeable with the standard glass door installed in the incubators you already have
- Recovery time drastically improved – The small access windows keep the changes in the environment inside the incubators to a minimum.

Temperature recovery data comparison of the regular glass door, optional shelf doors, and access windows

Test Conditions:

Room temperature: 25°C; Incubator temperature setting: 37.0°C; No gas control; The standard glass door, the bottom access window/shelf door left open for 30 seconds (temperature measured at 5 points on the upper shelf)

