



Incubate the Future with You

# Heat-Wise RH-3100

## “Where Innovation Meets Temperature”

Heat-Wise (RH-3100) with optional Hand Held Chamber minimizes the effects of environmental fluctuations.

Its straight forward and sturdy design guarantees user-friendliness without any complicated procedures.

Heat-Wise designed exclusively for the Cube(AD-3100) hand held chamber



## Safe, Simple and Robust

## Safe, Simple and Robust

Just position the hand-held chamber on the Heat-Wise and use it as a straightforward satellite incubator in your workstation.

Heat-Wise is offered for both the Split Lid Hand Held Chamber (HC3100-2) and the Single Lid Hand Held Chamber (HC-3100).

\*Compatible with existing AD-3100 models and no modification is required.

Mixed Gas to supply!  
RH-GP Optional  
Gas Port Adaptor




The back side



Lower: RH-3100  
Heat-Wise

Upper: HC-3100-2  
Hand Held Chamber (Split Lid)



## 【 Specifications 】

### Heat-Wise

Model	RH-3100
Size(mm)	W252×D270×H77mm
Temp. Control	Digital PID Control
Temp. control range	Room temp. +7°C~40°C
Temp. fluctuation range	±0.3
Temp. distribution accuracy	Aluminum Surface ±0.5°C
Gas connection port	IN/OUT (OD: 6mm)
Weight	3.2kg
Power	AC100-240V 1A(Max)50/60Hz

### Hand Held Chamber (Single Lid Design)

Model	HC-3100
Size(mm)	W180×D217×H20mm
Capacity	781mL
Dish number	35mm: 12 dishes 60mm: 6 dishes 4well: 4 dishes

### Hand-Held Chamber ( Split Lid Design)\*4 compartments

Model	HC-3100-2
One Compartment Size(mm)	W94×D100×H20mm
One Compartment Capacity	186 mL
Dish number for one compartment	35mm: 4 dishes 60mm: 1 dish 4well: 1 dish

